NTI-Kahla GmbH • Rotary Dental Instruments
Im Camisch 3 • D-07768 Kahla/Germany
Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
e-mail: nti@nti.de • www.nti.de

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GB Instructions for use CrownCut Sintered-diamond crown cutter

Read the instructions for use carefully before using the CrownCut.

The manufacturer does not accept any liability for damage caused by not adhering to the instructions for use or incorrect use of the product.

Contents

1. Indications of the CrownCut
2. Contraindication of the CrownCut
3. Recommended technique
4. Disinfection, cleaning and sterilization

1. Indications of the CrownCut:

The CrownCut can be used for cutting crowns fabricated from any type of material (porcelain, metals, zirconia).

Contraindication of the CrownCut

The CrownCut must not be used to prepare, separate or trim natural teeth.

3. Recommended technique:

Insert the crown cutter into the contra-angle of the micromotor.

The contra-angle should deliver a minimum of 50 ml of coolant water per minute. If necessary check this before use. Provide protection against water spray.

Important: It is essential to cut a deep "V" into the crown (see Fig. above). The cut is wider at the surface than deeper into the crown. The disc can become stuck between the two sections of the crown, bend and/or break if the patient moves. Not suitable for the lifting of crowns (risk of fracture).

The crown cutter should be cleaned with the cleaning stone G9920 after each use. If it is not cleaned, the diamond loses its cutting capacity. Adhere to the instructions for use supplied with the cleaning stone.

The crown cutting disc can be used until it is too small to be able to cut through the crown.

4. Disinfection, cleaning and sterilization:

Before using the CrownCut for the first time and immediately after every use disinfect it in a diamond instrument disinfectant (certified by VAH or EPA) according to the disinfectant manufacturer's instructions.

The CrownCut is sterilized in an autoclave at 134°C/5 min. or according to the unit manufacturer's instructions. Comply with DIN EN 13060 and DIN EN 285 (Packaging complies with EN ISO 11607/DIN EN 868-2ff.). **Observe guidelines from Robert Koch Institut (Hygiene Requirements in Dentistry) or the valid legal and hygiene regulations in the respective country!**

Immediately remove any instruments that are damaged, bent or no longer run concentrically. These instructions for use are available on our website at www.nti.de in the following languages: NL Dutch, GR Greek, HU Hungarian, SK Slovak, FI Finnish, PL Polish, PT Portuguese, SI Slovenian, CZ Czech, DK Danish, ER Estonian, LV Latvian, LT Lithuanian, SE Swedish.

Detailed recommendations for preparation in accordance with DIN EN ISO 17664

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