


Diamantinstrumente • Diamond Instruments • Instruments diamantés

801 (001)



Rund
Round
Rond




Größe/Size/taille ø 1/10 mm		007	008	009	010	012	014	016	018	021	023	025	029	035	042		
Länge-Length in mm - longueur																	
Shank	ISO	REF															
FGM	806 313 001 524...	801- ... M- FGM				012	014	018									
FGM	806 313 001 534...	801- ... C- FGM										023					
FG	806 314 001 524...	801- ... M- FG		007	008	009	010	012	014	016	018	021	023	025	029	035	042
FG	806 314 001 544...	801- ... SC- FG					012	014	016	018	021	023	025	029	035		
FG	806 314 001 534...	801- ... C- FG					012	014	016	018	021	023	025	029	035	042	
FG	806 314 001 514...	801- ... F- FG		007	008	009	010	012	014	016	018	021	023	025	029	035	
FG	806 314 001 504...	801- ... SF- FG					012	014	016	018	021	023	025	029	035		
FGXL	806 316 001 524...	801- ... M- FGXL				010	012	014	016	018	023						
FGXL	806 316 001 544...	801- ... SC- FGXL					014	016	018	023			029				
FGXL	806 316 001 534...	801- ... C- FGXL					014	016	018	023							
RA	806 204 001 524...	801- ... M- RA				010	012	014	016	018	021	023	029				

801L (697)



Rund, lang
Long Round
Rond, long



Größe/Size/taille ø 1/10 mm		012	014	016	023	029	035		
Länge-Length in mm - longueur									
Shank	ISO	REF							
FG	806 314 697 524...	801L- ... M- FG		012	014	016	023	029	035
FG	806 314 697 544...	801L- ... SC- FG			014	016	023	029	
FG	806 314 697 534...	801L- ... C- FG		012	014	016	023	029	035




801-023C-FGXL

802 (002)



Rund, mit Ansatz
Round with Collar
Rond avec extension



Größe/Size/taille ø 1/10 mm		008	009	010	012	014	016	018	023		
Länge-Length in mm - longueur		2,4	2,4	2,4	2,4	2,4	2,4	2,4	3,2		
Shank	ISO	REF									
FG	806 314 002 524...	802- ... M- FG		008	009	010	012	014	016	018	023
FG	806 314 002 544...	802- ... SC- FG			014	016	018	023			
FG	806 314 002 534...	802- ... C- FG				012	014	016	018	023	
FG	806 314 002 514...	802- ... F- FG				014	016	018			



801-023SC-FGXL

805 (010)



Umgekehrter Kegel
Inverted Cone
Cône renversé



Größe/Size/taille ø 1/10 mm		008	009	010	012	014	016	018	021	023		
Länge-Length in mm - longueur		1,0	1,0	1,0	1,2	1,6	1,6	1,6	1,9	2,5		
Shank	ISO	REF										
FGM	806 313 010 524...	805- ... M- FGM				012	014					
FG	806 314 010 524...	805- ... M- FG		008	009	010	012	014	016	018	021	023
FG	806 314 010 544...	805- ... SC- FG					014	016	018	021		
FG	806 314 010 534...	805- ... C- FG				012	014	016	018	021	023	
FG	806 314 010 514...	805- ... F- FG				014	016					
RA	806 204 010 524...	805- ... M- RA			010	012	014	016				

Bitte Drehzahlbereiche beachten (Drehzahltafel - Seite 11.12) / Please observe Recommended speeds (speed chart - page 11.12) / Respecter les vitesses de rotation recommandées (voir à la page 11.12)

M Ohne Ring/Without ring/Sans bague	Mittel/Medium/Standard	107-126 µm	ISO 524	F Roter Ring/Red ring/Bague rouge	Fein/Fine/Fin	40 µm	ISO 514
■ SC Schwarzer Ring/Black ring/Bague noire	Supergrob/Supercourse/Supergras	181 µm	ISO 544	■ SF Gelber Ring/Yellow ring/Bague jaune	Superfein/Superfine/Superfin	20 µm	ISO 504
■ C Grüner Ring/Green ring/Bague verte	Grob/Coarse/Gros	151 µm	ISO 534	□ UF Weißer Ring/White ring/Bague blanche	Ultrafein/Ultrafine/Ultrafin	15 µm	ISO 494