

Diamantstreifen mit und ohne Sägeverzahnung, Belegung einseitig Diamond finishing strips with and without serrated edge, single sided Strips diamantés avec ou sans dents de scie, couverte sur une seule face



diamantfreier Bereich / diamond free area
Partie sans diamantation (mm)

← 40 mm →



FS3-MX, FS3-FX, FS3-SFX

Farbkodierung/Colour Code/Code couleur

	Blue	Red	Yellow
Stärke/Thickness/grosueur (mm)	0,13	0,10	0,08
Breite/Width/largeur (mm)	3,75	3,75	3,75
Länge/Length/longueur (mm)	147,0	147,0	147,0



REF



FS3-MX




FS3-FX



FS3-SFX



FS3-X-SO 2 + 4 + 4



FS2-M, FS2-F, FS2-SF

← 16 mm →

Farbkodierung/Colour Code/Code couleur

	Blue	Red	Yellow
Stärke/Thickness/grosueur (mm)	0,13	0,10	0,08
Breite/Width/largeur (mm)	2,50	2,50	2,50
Länge/Length/longueur (mm)	147,0	147,0	147,0



REF



FS2-M



FS2-F



FS2-SF



FS3-M, FS3-F, FS3-SF

← 16 mm →

Farbkodierung/Colour Code/Code couleur

	Blue	Red	Yellow
Stärke/Thickness/grosueur (mm)	0,13	0,10	0,08
Breite/Width/largeur (mm)	3,75	3,75	3,75
Länge/Length/longueur (mm)	147,0	147,0	147,0



REF



FS3-M



FS3-F



FS3-SF